

Jobst UlcerCare Sizing Chart

Size	Ankle Circumferenc	Calf Circumference
Small	7" - 8 1/4" (17.8-21cm)	11 1/2" - 13 3/4" (29.2-34.9cm)
Medium	8 1/4" - 9 1/2" (21-24.1cm)	13 1/4" - 16" (33.7-40.6cm)
Large	9 1/2" - 10 3/4" (24.1-27.3cm)	15 3/4" - 18 1/2" (40-47cm)
X-L	10 3/4" - 12" (27.3-30.5cm)	18" - 21 1/4" (45.7-54cm)
XX-L	12" - 13" (30.5-33cm)	19" - 22" (48.3-55.9cm)
XXX-L	13" - 14" (33-35.6cm)	19 3/4" - 23" (50.2-58.4cm)
XXXX-L	14" - 15 1/4" (35.6-38.7cm)	20 1/2" - 23 3/4" (52.1-60.3cm)