

6856

Bariatric Waist-High Portable Medical Scale

Designed for obese and unsteady patients, this 1,000-lb capacity scale features a balanced, lightweight design. It's easily tipped back for transportation on its hard, large diameter, rubber wheels yet narrow enough to fit through standard doors. The battery-powered, LCD display features easy-to-use Body Mass Index calculation while weighing. The low-profile platform, convenient handrail, sturdy steel construction, and skid-proof mat add to patient assurance and stability.



SPECIFICATIONS

Model:	6856
Capacity:	1,000 lb x 0.2 lb / 450 kg x 0.1 kg
Platform Size:	24"W x 24"D / 61 cm x 61 cm
Display:	<ul style="list-style-type: none">• 0.75" high digits, high-contrast LCD• Body Mass Index calculation• LOCK button holds weight after weighing
Power:	6 "AA" size batteries (not included) or optional AC adapter.
Shipping Weight:	85 lb / 39 kg

Options include height rod and AC adapter.



Large LCD and Body Mass Index calculation