


Model: SB2

# TubBuddy

The TubBuddy is an ergonomic bathing transfer system designed to allow a patient to safely transfer over a bath tub without committing to any costly construction. This model is also capable of rolling over standard/elongated toilets to be used as a commode. The TubBuddy eliminates unnecessary transfers between toileting and bathing making it comfortable for both patient and caregiver meanwhile reducing costs and contributing to good patient hygiene.

## Additional product features include:

- 440 lb. weight capacity
- Rolls over standard/elongated toilets
- Tool-less assembly
- Low maintenance
- Removable / foldaway arms
- Safety mechanisms
- Corrosion, mold resistant
- Footrest
- 5" Locking castors
- Adjustable neck rest
- Lap / chest belts
- Commode system



- » ROLLING BASE 20"W X 28"D
- » BACK WIDTH 18"
- » BACK HEIGHT 19"
- » FLOOR TO SEAT HEIGHT 26"
- » FLOOR TO NECK REST HEIGHT 50"
- » MAX TOILET CLEARANCE 18"
- » BRIDGE LENGTH (STANDARD) 16"
- » INSIDE ARM TO ARM WIDTH 18"
- » OUTSIDE ARM TO ARM WIDTH 21"
- » SEAT DIMENSIONS 18"W X 17"D
- » TUB BASE (AT TOP) 19" W X 21" H
- » TUB BASE (AT BOTTOM) 16" W X 21" H
- » SEAT HOLE 7"W X 14"L
- » WEIGHT CAPACITY 440

