

Accurate. Quick. Easy to use.

Jamar[®] Plus+ Digital Goniometer.

The perfect goniometer for therapists who are looking for a ROM tool that is easier to use and more accurate than a standard plastic goniometer.

- Instant digital read-out to a tenth of a degree.
- "Zero" button allows you to zero the goniometer.

Start your measurement from any position and get an accurate measurement without needing to subtract the beginning value.

- Display numbers can be flipped with the press of a button, so you don't have to read upside-down measurements.
- "Hold" button freezes display for easier recording of results.
- Auto shut-off for longer battery life.
- Available in an 8" and 12½" versions.

Requires a CR2032 battery (included). One-year limited warranty (battery not covered).

0816-60-133 8" Goniometer \$34.95

0816-60-141 12½" Goniometer \$37.95



Patterson Medical
28100 Torch Parkway, Suite 700
Warrenville, IL 60555-3938

www.pattersonmedical.com

Call toll-free:
1-800-323-5547

Fax toll-free:
1-800-547-4333