



TILT TABLE

 350 lb capacity

Motorized, Trendelenberg

Item No.

TG2740.TB

FEATURES:

- 350 lb Weight Capacity
- Trendelenberg position 0 to -12 degrees
- UL Recognized Components
- Casters - Heavy Duty 5" Swivel with Locking Brakes
- Removable Footboard
- Safety Belts
- Angle Indicator
- Accommodates Patient Lifts
- 2" Upholstered Top
- Steel Frame

Width:	28"
Length:	78"
Height:	33"
Tilt Range:	0 to 85 degrees
Trendelenberg:	0 to -12 degrees
Base:	3" Tubular Steel
Finish:	Pearl Grey Standard; Porcelain Optional.
Upholstery:	Imperial Blue, Regimental Blue, Wedgewood, Blue Ridge, Dove, Graphite, Black, Burgundy, Mauve, Forest, Grotto, Sand, Camel.
Electrical:	120V, 6.5A, 780W, hospital grade plug.
Warranty:	Five (5) Year Limited Warranty.

*Nominal Dimensions Stated