

A Therapist
Favorite!



ComfySplints
We Make The World Comfy

Comfy™ Hand Thumb Orthosis

OUR HAND WRIST ORTHOSIS WITH EXTRA THUMB SUPPORT

SKIN SAFE

Anti-microbial and
moisture wicking
material for skin
breakdown
prevention



COMFORT

Super Soft and
removable
Headliner or
Terrycloth covers
are machine
washable

THUMB SUPPORT

Unique folding
thumb tabs allow
for extra thumb
control and support.



Shown in Navy Comfyrene version

ADJUSTABLE

Patented spine
bends to wrist
and forearm
needs

POSSIBILITIES

Comes in many variations:

- Comfyrene- Slim line, no cover
- Long Hand
- Long Arm
- Pediatric
- Broadcloth or Terrycloth cover
- 4 straps

RELAXING

An ideal
resting splint

SIDE WINGS

Prevents ulnar
and radial
deviation



Ambidextrous

	Wrist to Tip of Middle Finger	Total Length
Adult LG	8" and Above	13 ¾"
Adult	7" – 8"	12 ½"
Adult Sm	6" – 7"	11"
Pedi Large	4 ½" – 6"	9 ½"
Pedi Med	3" – 4 ½"	8"
Pedi Small	2 ½" – 3 ½"	7"

550990204	Comfy Hand Thumb
081621747	Comfyrene Hand Thumb
081566819	Comfy Hand Thumb- Broadcloth
081566801	Comfy HT -Broadcloth, Ad. Sm.

PATTERSON
MEDICAL

www.pattersonmedical.com

Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada Fax 1-877-566-8611