

# 3pp™ ThumSpica™ Plus



Wrist and thumb support includes a thermoplastic insert for added control.  
(Insert included for custom molding by a health care professional)



3pp™ ThumSling™



3pp™ ThumSling™ Long



3pp™ ThumSpica™



3pp™ Wrist Wrap™




3pp™ Elbow Wrap™



3pp™ Carpal Lift™



## 3pp™ ThumSpica™ Plus




**Fits circumference from**  
6" - 9" (15.2 - 22.8 cm)

### To Order

	<b>Right</b>	<b>Left</b>
One size	P3019-R	P3019-L

**Order Toll Free 888.378.7763** (U.S. only)  
Local/Int. 410.604.6393  
24 hour Fax 410.604.6398  
Mon. - Thurs. 8:30 - 5:30 pm, Fri. 8:30 - 5:00 pm EST



Remove insert prior to washing. Hand or machine wash soft wrap in soap and cool water. Wash insert in cool water and mild soap. Air dry. Frequent washing will extend product life.

Made in Mexico. Insert made in U.S.A.

Authorized European Representative  
 Emergo Europe  
Molenstraat 15  
2513 BH The Hague  
The Netherlands

©2009, 3-Point Products, Inc.

### Material

Nylon and spandex. Latex free  
Insert: Aquaplast low temperature thermoplastic.

[www.3pointproducts.com](http://www.3pointproducts.com)


# 3pp™ ThumSpica™ Plus



## Right Hand




1




2




3



4




5




6



7



8



9

Watch for our 3pp videos at [www.youtube.com/3pointproducts](http://www.youtube.com/3pointproducts)

3-Point Products, Inc

Stevensville, MD U.S.A

888-378-7763


# 3pp™ ThumSpica™ Plus



## Left Hand



1




2



3



4



5



6



7



8



9

Watch for our 3pp videos at [www.youtube.com/3pointproducts](http://www.youtube.com/3pointproducts)

3-Point Products, Inc

Stevensville, MD U.S.A

888-378-7763