



# Our Exclusive Pelvic Health Table

The first treatment table designed specifically for Pelvic Health Rehab Therapists,  
**it solves problems posed by traditional tables,**  
 namely uncomfortable positioning for the therapist and patient.

The Hi-Lo Table provides better internal evaluation and treatment of pelvic floor disorder (PFD), orthopedic conditions and bladder control—approximately 33 million Americans experience bladder problems each year. Other applications include gastrointestinal problems, sacroiliac dysfunction, infertility and coccyx pain. This innovative table allows therapists to achieve an ergonomic posture that greatly reduces stress on shoulders and spine.



**Adjustable Legrest**

**Adjustable Legrest** helps stabilize the patient, accommodating the patients height. Adjustments facilitate hip legs during an internal exam, and use the unsupported leg to activate the hip rotator.



**Side Cutouts**

**Side Cutouts** allow the clinician to get close to the patient, yet leave enough room for the patient to be comfortable. Standard treatment tables are wider, more difficult to get around. They can cause potentially damaging side-bend and rotation for therapists.



**Powered Center Section**

**Powered Center Section** lets the therapist effortlessly position and reposition the patient during treatment. Traditional tables often require therapists to use multiple pillows—a position that can cause back pain for the patient.



**Bar Activated**

**Bar-Activated**, Hi-Lo control enables clinicians to easily adjust table height anywhere around the table without interrupting treatment.

## Specifications:

Capacity	400 lbs.
Table Top	Five-piece. Overall 27"W, 72"L.
Head Section	15" x 13". Raises 25°. Lowers 45°. Contoured face/nose opening.
Armrests	Vertical 6" adjustments remain on a plane with the center section.
Center Section	19" x 27". Motorized to raise 20°.
Foot Section	38" x 27". Raises 75°. Dual gas spring supports.
Height Range	19"–38"

Call toll free  
 1-800-323-5547

In Canada, call 1-800-665-9200

Fax toll free  
 1-800-547-4333

In Canada, fax 1-800-566-8611

  
**PATTERSON  
 MEDICAL**  
 A Patterson Medical Exclusive