

Bariatric Patient Lifter Rhino P500



Heavy Duty
Hydraulic Cylinder

- Deluxe patient lifter with powder coated carbon steel constructed frame with high stress areas reinforced.
- Heavy duty hydraulic cylinder.
- Adjustable base.
- 3" casters with parking locks.
- Weight limit: 1000 lbs.

MODEL	SLING HEIGHT	BASE CLEARANCE	BASE WIDTH	BASE LENGTH	OVERALL HEIGHT	WEIGHT W/O SLING	PACKAGING
P500	0" - 42"	5"	57"	58"	84" - 91"	150 lbs.	1 / CASE