

# Jogger

2 year warranty

## Supportive Style On-the-Go

- Designed for parents and families who want to take long strolls through the city or through rough trails
- Easy to fold and transport
- Lightweight aluminum frame
- Includes 8-point trunk / pelvic harness system

Special Tomato® Jogger - 62000000	
Stroller Only Seat to Foot Well	14"
Liners Size 1 and 2 Seat to Foot Well	16"
Sitter Size 1 Seat to Foot Well	17"
Sitter Size 2 Seat to Foot Well	18"
Open Chair Dimensions	43"H × 26"W × 42"D
Folded Chair Dimensions	16"H × 33"W × 21"D
Handle Height to Floor	35" - 43"
Seat	14"W × 12"D
Back Height	25"
Recline Range	100° - 140°
Max. User Weight	110 lbs.
Stroller Weight	27 lbs.



# Features & Accessories



- 1 Size 1 and 2 Soft-Touch® Liners**  
Attach to the stroller, with color choice of Cherry, Lilac, Aqua or Dark Gray. Sold separately, see page 16.
- 2 Size 1 and 2 Soft-Touch® Sitters**  
Attach to the stroller, with color choice of Lilac, Aqua or Dark Gray. Sold separately, see page 9.
- 3 Added Comfort**  
Padded upholstery with reclining back support, functional rear suspension and foot well with side guard.
- 4 Pneumatic Wheels**  
Front 12" pneumatic swivel wheel with lock-out and quick-release 12" rear pneumatic wheels.
- 5 One-step, compact folding design**  
Canopy folds with Jogger.
- 6 Sun Shade**  
Removable canopy with window and height adjustable handle.
- 7 Medical Necessities Bag**  
62010000  
Sold separately.
- 8 Head Support**  
62020000  
Sold separately.
- 9 Rain Cover**  
62030000  
Sold separately.