

Dorsal Carpal Tunnel Splint

US Patent #6,773,410

by **RCAI**



44DCT



Indicated for carpal tunnel syndrome, ten-donitis, lateral epicondylitis, wrist sprain or fracture, the Dorsal Carpal Tunnel Splint limits wrist flexion and extension during repetitive hand motion. Made of lightweight, yet durable Kydex® plastic, the splint can be heat molded for patient customization and positioning. Hand wash or gentle machine wash liner. Tumble dry.

- Lined with breathable Orthowick™ to absorb perspiration and moisture.
- Narrow palmar roll supports the hand while fingers remain free for grasp.
- Prevents damaging wrist flexion, restoring finger sensation and hand strength.
- Wear while typing, writing, or any repetitive work related activity.
- Comfortable enough to be worn at night during rest.

Size	Measure width of palm	Model # Right	Model # Left
Ex Small	2" – 2 1/2" (5.1 – 6.4cm)	44DCT-ES-R	44DCT-ES-L
Small	2 1/2" – 3" (6.4 – 7.6cm)	44DCT-S-R	44DCT-S-L
Medium	3" – 3 1/2" (7.6 – 8.9cm)	44DCT-M-R	44DCT-M-L
Large	3 1/2" – 4" (8.9 – 10.2cm)	44DCT-L-R	44DCT-L-L
X-Large	4" – 4 1/2" (10.2 – 11.4cm)	44DCT-XL-R	44DCT-XL-L




Restorative Care of America
INCORPORATED

Phone: (727) 573-1595 or (800) 627-1595
Fax: (727) 573-1886 or (800) 545-7938
12221 33rd Street N, St Petersburg, FL 33716
www.rcai.com

Hand-Crafted and Made In The USA

Manufacturer of the Original
and ONLY Multi Podus® System



Limited Warranty Coverage: Any defect in material or workmanship solely limited to the Kydex® plastic frame or metal portion or hinge mechanisms for 1 (one) year from date of purchase.

Dorsal Carpal Tunnel Splint

US Patent #6,773,410

by **RCAI**



44DCT



Indicated for carpal tunnel syndrome, ten-donitis, lateral epicondylitis, wrist sprain or fracture, the Dorsal Carpal Tunnel Splint limits wrist flexion and extension during repetitive hand motion. Made of lightweight, yet durable Kydex® plastic, the splint can be heat molded for patient customization and positioning. Hand wash or gentle machine wash liner. Tumble dry.

- Lined with breathable Orthowick™ to absorb perspiration and moisture.
- Narrow palmar roll supports the hand while fingers remain free for grasp.
- Prevents damaging wrist flexion, restoring finger sensation and hand strength.
- Wear while typing, writing, or any repetitive work related activity.
- Comfortable enough to be worn at night during rest.

Size	Measure width of palm	Model # Right	Model # Left
Ex Small	2" – 2 1/2" (5.1 – 6.4cm)	44DCT-ES-R	44DCT-ES-L
Small	2 1/2" – 3" (6.4 – 7.6cm)	44DCT-S-R	44DCT-S-L
Medium	3" – 3 1/2" (7.6 – 8.9cm)	44DCT-M-R	44DCT-M-L
Large	3 1/2" – 4" (8.9 – 10.2cm)	44DCT-L-R	44DCT-L-L
X-Large	4" – 4 1/2" (10.2 – 11.4cm)	44DCT-XL-R	44DCT-XL-L




Restorative Care of America
INCORPORATED

Phone: (727) 573-1595 or (800) 627-1595
Fax: (727) 573-1886 or (800) 545-7938
12221 33rd Street N, St Petersburg, FL 33716
www.rcai.com

Hand-Crafted and Made In The USA

Manufacturer of the Original
and ONLY Multi Podus® System



Limited Warranty Coverage: Any defect in material or workmanship solely limited to the Kydex® plastic frame or metal portion or hinge mechanisms for 1 (one) year from date of purchase.