


For below knee make selection by measuring ankle, and calf circumference, and leg length from heel at floor to bend in knee. For thigh high make selection by measuring ankle, calf, and thigh circumference, and leg length from heel at floor to mid thigh.

20-30 and 30-40 Stockings Size Chart

Size	circumference			length
	ankle	calf	thigh	to thigh
Small	7" - 8 1/4"	11" - 14"	up to 20"	up to 26"
Medium	8 1/8" - 9 3/8"	13 1/2" - 16"	up to 22"	up to 28"
Large	9 1/4" - 11"	15 1/2" - 18"	up to 23"	up to 29"
X-Large	11 1/8" - 12 1/8"	17 1/2" - 20"	up to 24"	up to 31"
2X-Large	12 1/2" - 13 3/4"	19 1/2" - 22"	N/A	N/A
3X-Large	12 1/2" - 13 3/4"	21 1/2" - 24"	N/A	N/A
	Regular: greater than 15"		Short length below knee stockings are only available in small through x-large sizes.	
	Short: 15" or less			