

The GuardION™ Advantage

Fabrifoam's® innovative technology infuses GuardION orthotics with silver and natural minerals, a compound of micro-crystals with the ability to reflect and emit beneficial anions and far infrared rays. When applied, body heat and motion activate the compound to release anions and FIR back to the skin, stimulating the body to boost blood flow, oxygenation and accelerate the natural healing process.

Dynamic compression and support with enhanced therapeutic properties to help reduce pain, swelling and inflammation, suppress bacterial growth, minimize odors and aid with moisture evaporation.

Made in the USA
Manufactured by
Fabrifoam® Products
900 Springdale Drive
Exton, PA 19341
www.Fabrifoam.com
800-577-1077



The Original Therapeutic Material®

CarpalGard™

Recommended for carpal tunnel syndrome (CTS), wrist pain, wrist strain, arthritis, tendinitis, UCL and other hand and wrist ailments. Moderate to semi-rigid support. Non-bulky and flexible for preventive care and rehab. Encourages range of motion. Adjustable compression and custom fit.



GUARDION™

Advanced Technology in Far Infrared & Anion Therapy

CarpalGard™

Your CarpalGard™ is made from Fabrifoam's® triple-layer X-treme™ material to provide the ultimate in support and cushioning with the power of a rigid brace and the comfort of soft-bracing. This fabric/foam composite is latex-free, neoprene-free, non-slip, breathable and moisture-wicking for safe and effective musculoskeletal care. Infused with GuardION™ for the added benefits that negative ions and FIR can deliver.

Washing Your Fabrifoam®

Close fastener(s). Hand rinse in warm water after use. Hand wash with a mild fragrance-free, dye-free soap as needed. Rinse well and gently squeeze out excess water. Hang or lay flat, foam side up, to air dry. Avoid use of bleach or fabric softener. Do not machine wash, machine dry or dry clean. This may damage your product.

If sensitivity or discomfort occurs, discontinue use and consult your healthcare professional.

Basic Application

Note: Photos show right-hand application; reverse for left hand.

Step 1:



Place thumb through slot (foam side to skin) with short strap positioned on the back of hand.

Step 2:



Wrap long strap snugly around and over short. Secure with internal hook & loop dot.

Step 3:



Continue wrapping strap snugly once around hand.

Step 4:



Continue wrapping around wrist.

Step 5:



Secure with hook & loop tab.

For comfortable support do not over tighten. Re-wrap to adjust compression if needed.