


Product Specification Sheet



Specifications:

Vendor Parts Number: M30, M50 & M70.

Description: 1:1 weight-to-resistance ratio, TIG welded high strength steel.

Starting resistance: 1.25 lbs.

Adjustability: 2" interval rope adjustment bar w/ spring-loaded plungers.

Split spacer, high end ball bearing wheels.

Standard: 1x pulley handle & 3x 0.5 lb. accessory weights.

The pulleys can be wall mounted, module mounted or free-standing.

Assembly: None required.

Product Dimensions: 84" (T) x 18" (W) x 8" (D). 70 lbs. / 90 lbs. / 110 lbs.

Shipping Dimensions: 87" (L) x 18" (W) x 14" (H), 150 lbs. (Class 85).

SHIPS VIA TRUCK 

Application: To restore functional qualities used in activities of daily living.

Functions:

The RehabPro Standard Pulleys, with their accessories, are considered one of the most versatile pieces of rehab equipment to improve tissue tolerance to physical load. Rope adjustability, low starting weights and flexible base options offer the greatest number of exercise in the least amount of space. Standard pulleys have restored function to thousands of clients on multiple continents for almost a half a century.