


Product Specification Sheet



Specifications:

Vendor Parts Number: M100S

Description: 6:1 weight-to-resistance ratio, TIG welded high strength steel. Resistance as low as 0.85 lbs.

Adjustability: 2" interval rope adjustment bar w/ spring-loaded plungers.

Split spacer, high end ball bearing wheels.

Standard: 2x pulley handle.

The pulleys can be wall mounted, module mounted or free-standing.

Assembly: None required.

Product Dimensions: 84" (T) x 18" (W) x 8" (D). 175 lbs.

Shipping Dimensions: 87" (L) x 18" (W) x 14" (H), 250 lbs. (Class 70).

SHIPS VIA TRUCK 

Application: To restore functional qualities used in activities of daily living.

Functions:

The RehabPro Speed Pulleys create additional exercise dosage adjustments not possible with standard pulleys. Adjustment of speed and application of multi-directional force take stimulation of soft tissue to a new level improving functional qualities needed in activities of daily living.