

Product Specification Sheet



(Standard Accessories: 2x Pulley Handles)

Functions:

This bilateral pulley widens the scope of rehab by offering the option of bilateral work. Bilateral work could mean doing isotonic work with two extremities or superimpose dynamic movement of one extremity on a resisted holding pattern of the spine. This pulley can behave as a standard pulley but at a 50% reduction in resistance and a reduced distance of rope travel. M100B can be wall or module mounted.

Specifications:

Vendor Parts Number: M100B

Construction:

2:1 weight-to-resistance ratio physical weights stack. The weight stack consists of 20x 5.0 lb. weights which translates into a resistance of 2.5 lbs. per weight plate.

The pulley is TIG welded and made of durable, powder coated high strength steel

The 3/16" non-stretchable marine rope adjusts in height and length along an 84.5" - 2" interval rope adjustment bar w/ spring-loaded plungers.

The pulleys comes with 3/4" solid steel, chromium plated, self-adjusting guide rods ensuring smooth tracking of the weight stack.

Assembly: None required.

Product Dimensions: 84" (T) x 18" (W) x 8" (D), 150 lbs.

Shipping Dimensions: 87" (L) x 18" (W) x 14" (H), 200 lbs. (Class 70)

SHIPS VIA TRUCK 

Warranty:

1 year bumper to bumper all parts.
5 years on all moving parts,
Life time on frame.