

TributeWrap™ **NEW**

Adjusts to your patient and their lifestyle.



TributeWrap™

Adjustable variable compression foam garments for lymphedema and edema.

TributeWrap™ is an off-the-shelf garment that comfortably adjusts to your patients' limb and lifestyle. This versatile garment can be used during the evening, night, or low activity — such as working at a desk job or resting on the couch.

While your patients rest, its proven foam technology works hard to break up hardened or fibrotic tissue to help support soft and healthy skin care. Most importantly, it allows your patients to maintain their gains made during therapy, while also allowing the flexibility needed to support further reduction.

- **Durable straps** help your patients adjust their garment for the right fit
- **Donning tabs** make it easy to get a better grip and pull for improved contour to the limb
- **Powered by Coolcore®** to rapidly wick away sweat and moisture and regulate body temperature
- **Chevron channeling** is designed based on MLD principles of directing fluid to collateral pathways
- **Low profile height** helps provide a comfortable garment experience

TributeWrap™ Glove

Size	Circumference (cm)		Orientation	Item No.
	Palm (B)	Calf (C)		
Small	16–20	16.5–21	Left	7101702
			Right	7101705
Medium	18–22	17.5–22	Left	7101703
			Right	7101706
Large	20–24.5	18.5–23	Left	7101704
			Right	7101707

TributeWrap™ Wrist to Axilla

Size	Circumference (cm)			Arm Length (C–G)	Orientation	Item No.
	Wrist (C)	Elbow (E)	Axilla (G)			
Small	16–21	22–28	25–31	43–47	Left	7101695
					Right	7101698
Medium	17–22	26–32	30–37		Left	7101696
					Right	7101699
Large	18–23	30–36	35–44		Left	7101697
					Right	7101700

TributeWrap™ Below Knee

Size	Circumference (cm)		Length (cm)		Orientation	Item No.
	Ankle (B)	Calf (D)	Foot (A–Y*)	Calf (Y*–D)		
Small	24–30	30–36	19–22	40–44	Left	7101709
					Right	7101712
Medium	27–33	33–41.5			Left	7101710
						Right
Large	30–36	36–45			Left	7101711
						Right

*Y measurement is at the base of the heel.

Effective, easy-to-use, and comfortable —
a great option for your patient!