



## TributeWrap™ **NEW**

Adjustable variable compression foam garments for lymphedema and edema.



# Explore TributeWrap.™ **NEW**

TributeWrap is an off-the-shelf adjustable foam garment for lymphedema or edema management during the evening, night or low activity. TributeWrap comfortably adjusts to your patients' unique shape and lifestyle to help maintain their gains made during therapy, while allowing the flexibility to support further reduction.




## Durable straps and easy-to-use donning tabs

Adjusts to your patient's shape or size with donning tabs that make it easy to get the right fit and contour.



## Angled straps and darting at flexion points

Allows freedom to move to any resting or sleeping position.



POWERED BY  **coolcore**

Innovative fabric with chemical-free Cooling Power that rapidly wicks away sweat and moisture to regulate body temperature.



## Chevron channeling

Designed based on MLD principles of directing lymphatic fluid to collateral pathways.



## Low profile height

Slim design benefits patient comfort.




## Proven foam technology

Helps soften and break up hardened or fibrotic tissue.

## Quick and easy fitting

Only 2–4 measurements are needed to select size for each garment.



## Sizing

### TributeWrap Glove

Size	Circumference	
	Palm	Wrist
Small	16–20 cm	16.5–21 cm
Medium	18–22 cm	17.5–22 cm
Large	20–24.5 cm	18.5–23 cm

### TributeWrap Wrist to Axilla

Size	Circumference			Arm Length (C–G)
	Wrist	Elbow	Axilla	
Small	16–21 cm	22–28 cm	25–31 cm	43–47 cm
Medium	17–22 cm	26–32 cm	30–37 cm	
Large	18–23 cm	30–36 cm	35–44 cm	

### TributeWrap Below Knee

Size	Circumference		Foot Length (A–G)	Calf Length (B–D)
	Ankle	Calf		
Small	24–30 cm	30–36 cm	19–22 cm	40–44 cm
Medium	27–33 cm	33–41.5 cm		
Large	30–36 cm	36–45 cm		



### TributeWrap™ Glove

Size	Style	Item No.
Small	Left	136783
	Right	136786
Medium	Left	136784
	Right	136787
Large	Left	136785
	Right	136788

### TributeWrap™ Wrist to Axilla

Size	Style	Item No.
Small	Left	136777
	Right	136780
Medium	Left	136778
	Right	136781
Large	Left	136779
	Right	136782

### TributeWrap™ Below Knee

Size	Style	Item No.
Small	Left	136789
	Right	136792
Medium	Left	136790
	Right	136793
Large	Left	136791
	Right	136794

For more information or to find a dealer near you, please contact:

**L&R USA INC.**  
 3880 W Wheelhouse Road  
 Milwaukee, WI 53208 USA

Phone: 855-892-4140, 414-892-4140  
 Fax: 414-892-4150  
 Email: [inquiries@us.LRmed.com](mailto:inquiries@us.LRmed.com)  
 Website: [www.Lohmann-Rauscher.us](http://www.Lohmann-Rauscher.us)

Coolcore® is a registered trademark of Coolcore, LLC.