

# Electrotherapy

## InTENSity™ Twin Stim® III

The InTENSity™ Twin Stim® III combine TENS and EMS therapy in one easy-to-use digital device. When operating in TENS mode the unit has four modes and three modes when operating in EMS mode. An industry leader at 105 mA, this device comes standard with an AC adaptor.



### Product Features:

- Simple to operate, portable device for pain management
- Unit combines TENS and EMS therapy for optimal treatment
- Four modes for TENS therapy: Burst, Normal, Mod 1 and Mod 2
- Three modes for EMS therapy: Synchronous, Asynchronous and Delay
- Comes complete with carrying case, two sets of lead wires, 9-volt battery, AC adaptor, manual and pack of four self-adhesive reusable electrodes
- Power Source: Standard 9-volt battery or AC adaptor
- Dimensions: 4.5" x 2.55" x 0.9"
- Weight: 0.28 lb. with battery
- One year warranty
- HCPCS Code: E0730\* for TENS Unit
- HCPCS Code: A4595\* for electrotherapy supplies



### Product Specifications:

	TENS	EMS
Wave Form	Monophasic Square	Biphasic
Pulse Rate	1-150 Hz	1-150 Hz
Pulse Width	50-300 $\mu$ s	50-300 $\mu$ s
Output Voltage	Constant Current	Constant Current
Current	0-105 mA	0-105 mA
Modes	4	3
Timer	1-60 Constant	1-60 Constant

Item #	Description	Qty
DI3717	InTENSity™ Twin Stim® III	each

\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**

(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com