

SUPPORT PRO 7 SERIES

HCPCS CODING: E2602
GENERAL USE




PRODUCT FEATURES

Kwik Strap loops around vertical canes to securely attach cushion to wheelchair during transfers.


Non-countoured seating surface allows for easy transfers.

Coccyx relief area unloads bony prominences.



Waterfall front eliminates pressure points in leg trough.

750 lb. weight capacity foam provides increased and reduces potential high-pressure areas



SEE CHART FOR DIMENSIONS

STANDARD COVER



COMFORT-TEK FABRIC: the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.



OPTIONAL QUADRAGEL PACK

The QuadraGel pack regulates body temperature and redistributes pressure for comfort.



SUPPORT PRO 7 SERIES

WITHOUT QUADRAGEL

| comfort-tek part # | w x l | front | rear | wt. limit |
|--------------------|-----------|-------|------|-----------|
| 53-2218-B | 22" x 18" | 3.5" | 3.5" | 750 lbs |
| 53-2220-B | 22" x 20" | 3.5" | 3.5" | 750 lbs |
| 53-2418-B | 24" x 18" | 3.5" | 3.5" | 750 lbs |
| 53-2420-B | 24" x 20" | 3.5" | 3.5" | 750 lbs |
| 53-2618-B | 26" x 18" | 3.5" | 3.5" | 750 lbs |
| 53-2620-B | 26" x 20" | 3.5" | 3.5" | 750 lbs |
| 53-2818-B | 28" x 18" | 3.5" | 3.5" | 750 lbs |
| 53-2820-B | 28" x 20" | 3.5" | 3.5" | 750 lbs |
| 53-3018-B | 30" x 18" | 3.5" | 3.5" | 750 lbs |
| 53-3020-B | 30" x 20" | 3.5" | 3.5" | 750 lbs |

NOTES

+ CUSHIONS AVAILABLE WITH RIGID INSERTS & EXTRA INCONTINENT PROOF LINER. SEE PAGE 19 FOR DETAILS

SUPPORT PRO 7 SERIES

WITH QUADRAGEL

| comfort-tek part # | w x l | front | rear | wt. limit |
|--------------------|-----------|-------|------|-----------|
| 53G-2218-B | 22" x 18" | 3.5" | 3.5" | 750 lbs |
| 53G-2220-B | 22" x 20" | 3.5" | 3.5" | 750 lbs |
| 53G-2418-B | 24" x 18" | 3.5" | 3.5" | 750 lbs |
| 53G-2420-B | 24" x 20" | 3.5" | 3.5" | 750 lbs |
| 53G-2618-B | 26" x 18" | 3.5" | 3.5" | 750 lbs |
| 53G-2620-B | 26" x 20" | 3.5" | 3.5" | 750 lbs |
| 53G-2818-B | 28" x 18" | 3.5" | 3.5" | 750 lbs |
| 53G-2820-B | 28" x 20" | 3.5" | 3.5" | 750 lbs |
| 53G-3018-B | 30" x 18" | 3.5" | 3.5" | 750 lbs |
| 53G-3020-B | 30" x 20" | 3.5" | 3.5" | 750 lbs |

NOTES

+ CUSHIONS AVAILABLE WITH RIGID INSERTS & EXTRA INCONTINENT PROOF LINER. SEE PAGE 19 FOR DETAILS