

SUPPORT PRO

TOP
SELLER



ZERO ELEVATION

The Support Pro Zero Elevation allows for greater freedom of movement for standard positioning issues.



ANTI-THRUST

The Support Pro Anti-Thrust provides anterior buildup at midline to prevent forward ischial migration.



WEDGE

The Support Pro Wedge controls forward leaning by adding slight anterior height.

ABOUT THE SUPPORT PRO

- Non-contoured seating surface allows for easy transfers.
- Optional QuadraGel pack insert.
- Waterfall front eliminates pressure points in leg trough.
- Large radius bottom provides versatile surface for sling or flat seated wheelchair.
- Kwik Strap loops around vertical canes to securely attach cushion to wheelchair during transfers.
- Coccyx relief area unloads bony prominences.
- Revolutionary Moleculon® foam provides increased comfort and reduces potential high-pressure areas.

QUADRAGEL® PACK

- The QuadraGel pack regulates body temperature and redistributes pressure for comfort.



WEIGHT CAPACITY

| | | | |
|----------|---------|----------|---------|
| 16" Wide | 250 lbs | 20" Wide | 350 lbs |
| 18" Wide | 300 lbs | 22" Wide | 400 lbs |

STANDARD SIZES

| | 16" Wide | 18" Wide | 20" Wide | 22" Wide |
|----------|----------|----------|----------|----------|
| 16" Deep | 16 x 16 | 18 x 16 | | |
| 18" Deep | 16 x 18 | 18 x 18 | 20 x 18 | 22 x 18 |

ZERO ELEVATION

| W x D | Zero Elevation | Zero Elevation w/Gel | Zero Elevation w/Abductor | Zero Elevation w/Abductor & Gel |
|-----------|----------------|----------------------|---------------------------|---------------------------------|
| 16" x 16" | 52A | 52GA | 52PA | 52PGA |
| 16" x 18" | 52A18 | 52GA18 | 52PA18 | 52PGA18 |
| 18" x 16" | 52B16 | 52GB16 | 52PB16 | 52PGB16 |
| 18" x 18" | 52B | 52GB | 52PB | 52PGB |
| 20" x 18" | 52C | 52GC | 52PC | 52PGC |
| 22" x 18" | 52D | 52GD | 52PD | 52PGD |

WEDGE

| W x D | Wedge | Wedge w/Gel | Wedge w/Abductor | Wedge w/Abductor & Gel |
|-----------|-------|-------------|------------------|------------------------|
| 16" x 16" | 55A | 55GA | 55PA | 55PGA |
| 16" x 18" | 55A18 | 55GA18 | 55PA18 | 55PGA18 |
| 18" x 16" | 55B16 | 55GB16 | 55PB16 | 55PGB16 |
| 18" x 18" | 55B | 55GB | 55PB | 55PGB |
| 20" x 18" | 55C | 55GC | 55PC | 55PGC |
| 22" x 18" | 55D | 55GD | 55PD | 55PGD |

ANTI-THRUST

| W x D | Anti-Thrust | Anti-Thrust w/Gel | Anti-Thrust w/Abductor | Anti-Thrust w/Abductor & Gel |
|-----------|-------------|-------------------|------------------------|------------------------------|
| 16" x 16" | 57A | 57GA | 57PA | 57PGA |
| 16" x 18" | 57A18 | 57GA18 | 57PA18 | 57PGA18 |
| 18" x 16" | 57B16 | 57GB16 | 57PB16 | 57PGB16 |
| 18" x 18" | 57B | 57GB | 57PB | 57PGB |
| 20" x 18" | 57C | 57GC | 57PC | 57PGC |
| 22" x 18" | 57D | 57GD | 57PD | 57PGD |