

Adjustable glideboard promotes optimal spinal alignment and training variations.

Easy Two-Step Folding and transport wheels for moving unit.

LAT Bars promote line-of-pull specificity for arm-cable exercises.

Center-pulley Attachment Positions offer resistance challenges at six levels.

Automatic Rail Locks allow for explosive plyometric jumping.

Hydraulic Rail Lift assists level changes and easy two-step folding.

PRODUCT SPECIFICATIONS

IN USE	108" x 38" x 69" [L/W/H] (2.7 m x 1.0 m x 1.8 m)
FOLDED	20" x 38" x 56" [L/W/H] (.5 m x 1.0 m x 1.4 m)
UNIT WEIGHT	180 lbs (81.65 kg)
RAIL CONSTRUCTION	Steel reinforced extruded aluminum rails
RESISTANCE LEVELS	1% to 62% of bodyweight 22 calibrated incline levels, with release lever to lock rails into place
SQUAT STAND	Three height levels
STORAGE	Hydraulic lift assists folding; Rolls upright when folded
USER HEIGHT	Up to 6'10" (2.1 m)
CAPACITY	Maximum total 650 lbs (295 kg)
WARRANTY	Frame: 5 years Parts & Upholstery: 1 year Foam & Rubber: 90 days

The Recovery Series Encompass™ provides the optimum functional training experience. Works multiple muscle groups together for a complete workout that challenges proprioception and core stabilizers by integrating multi-plane movement and unrestricted range of motion. Perfect for your physical therapy wellness program.

1. Now included with the unit are Closed Chain Platform and BAPS Adapter, Retractable Dip Bars, Squat Handle Bars, Slide Distance Regulator.
2. Adjustable glideboard promotes optimal spinal alignment and training variations.
3. Automatic rail locks allow for explosive plyometric jumping.
4. Hydraulic rail lift assists level changes and easy two-step folding.
5. Multiple center-pulley attachment positions double arm-cable load and increases lower extremity ROM.
6. LAT Bars promote line-of-pull precision for arm-cable exercises.
7. Clinical Complete Accessory Package available for purchase. Package includes: Leg Pulley System, Press Bar, SCRUNCH® and Weight Bar.
8. Additional accessories available for purchase include: Pull-Up Bar (3Grip), Standing Platform and Telescoping Toe Bar.