



ProBasics K7


HEAVY DUTY WHEELCHAIR



The **ProBasics K7 Heavy Duty Wheelchair** features a robust design to accommodate patients with wider seat and higher weight capacity requirements. With seat widths from 24 to 28 inches, a weight capacity of 450 - 600 pounds and a heavy-duty slide tube frame, the ProBasics K7 offers a rugged yet attractive option for the bariatric patient population.



Removable, desk length, padded arms



Heavy duty vinyl upholstery



Aluminum footplate

FEATURES

- Heavy Duty, slide tube silver vein steel frame
- 8" x 2" semi-pneumatic casters are adjustable in 3 positions
- Dual axle provides easy transition of seat height to hemi-level
- 24" rear urethane wheels with durable hand rims
- Removable desk length arms with padded armrests
- Heavy duty vinyl padded upholstery with chart pocket
- Tool-free footrests include heel loops and Aluminum footplates
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite elevating positions and aluminum footplates
- Easy push-to-lock wheel locks
- HCPCS* Code: K0007
K00195 or E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	22" 24" 26" 28"
Seat Depth:	18" 20"
Seat-to-Floor Height ¹ :	Adult 21" Hemi 19"
Back Height:	16" 18"
Back Style:	Fixed
Overall Width w/Arms:	30" 32" 34½" 37"
Overall Depth:	33" 34" (without riggings)
Weight Capacity:	450 - 600 lb
Product Weight ² :	51 - 64 lb
Carton Dimensions:	35" x 13" x 38"
Shipping Weight:	58 - 70 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	Weight Limit	UPC	UOM
WC72218DS	22" x 18"	Removable, Desk-Length	Swing-Away Footrests	450 lb	815067070855	each
WC72418DS	24" x 18"	Removable, Desk-Length	Swing-Away Footrests	450 lb	815067070862	each
WC72620DS	26" x 20"	Removable, Desk-Length	Swing-Away Footrests	500 lb	815067072057	each
WC72820DS	28" x 20"	Removable, Desk-Length	Swing-Away Footrests	600 lb	815067072064	each
WC72218DE	22" x 18"	Removable, Desk-Length	Elevating Legrests	450 lb	815067070879	each
WC72418DE	24" x 18"	Removable, Desk-Length	Elevating Legrests	450 lb	815067070886	each
WC72620DE	26" x 20"	Removable, Desk-Length	Elevating Legrests	500 lb	815067072071	each
WC72820DE	28" x 20"	Removable, Desk-Length	Elevating Legrests	600 lb	815067072088	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.