


800-525-9339 • 651-674-8301
www.swedeo.com

- Light support for mild injuries, aching and/or weakened elbows
- Seamless design for a comfortable fit
- Gentle compression helps reduce inflammation
- Soft, ventilated knit material is suitable for extended use
- Contours to the natural shape of the elbow for comfort
- Bilateral design fits right or left elbow

MODEL #	SIZE	MEASUREMENT	
ELB-6517-SML	SM	7 3/4-9"	20-23cm
ELB-6517-MED	MED	9-10 1/4"	23-26cm
ELB-6517-LRG	LRG	10 1/4-11 1/2"	26-29cm
ELB-6517-1XL	XL	11 1/2-12 1/2"	29-32cm
Your Elastic Support should be a snug and comfortable fit. Not too tight, not too loose. If in doubt, try next larger size.			



Measure around elbow joint with arm extended

Functional & Comfortable

The Swede-O™ Elastic Elbow Sleeve provides light support for mild injuries, aching and/or weakened elbows. Gentle compression helps reduce inflammation. The seamless, ventilated material contours to the natural shape of the elbow for comfort during extended use. Discreet under clothing; bilateral design fits right or left elbow.



Specifications

Product Selection:

ELB-6517 Elastic Elbow Sleeve

Product Weight: .3 lb

Dimensions: 5.5" x 4.5" x 1.5"

Case Weight: 5 lb

Case Dimensions: 12" x 10" x 8"

Case Quantity: 12

Materials: Polyester, Rubber, Nylon

Available From:

