

AchilloTrain®

ACTIVE SUPPORT FOR ACUTE RELIEF OF THE ACHILLES TENDON.



viscoelastic insert



INDICATIONS

- Achilles tendinitis
- Achilles tendon partial ruptures
- Post-procedure (complete tendon ruptures)
- Swelling and pain
- Bursitis
- Haglund's deformity
- Need for tension reduction of the Achilles tendon (heel cup)

FEATURES & BENEFITS

- Three-dimensional Train active knit provides medical grade therapeutic compression. Anatomically contoured design promotes all-day wearing comfort.
- Viscoelastic inserts combined with medical grade compression support lateral and medial stabilization, control edema and reduce pain.
- Removable heel cup relieves the Achilles tendon.
- Separate heel cup is provided for unaffected leg to offset the length difference.
- Designed to support ankle at a 90° angle for acute relief (recent injury or surgery). Ideal for sitting and walking.
- Medical grade compression promotes increased proprioception and joint stabilization.
- Right and left specific for optimum fit.



Size	1	2	3	4	5	6	Order number
Circumference in inches	6 3/4 - 7 1/2	7 1/2 - 8 1/4	8 1/4 - 9	9 - 9 7/8	9 7/8 - 10 5/8	10 5/8 - 11 3/8	110110130*0+0°
Circumference in cm	17-19	19-21	21-23	23-25	25-27	27-29	

*1 = nature, *7 = black, *8 = titan, +6 = right, +7 = left, ° = size

