

Bariatric Transfer Bench



FIGURE A



FIGURE B

- Tool-free assembly of back, legs and arm
- A-frame provides additional stability
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- All transfer benches are reversible to accommodate any bathroom
- Euro-style clip ensures a secure fit (Figure A)
- No exposed hardware to injure patient
- Large, locking suction cups for added safety (Figure B)
- Weight Capacity 500 lbs

PRODUCTS

Description

12025KD-2	Knocked Down, Tool-free Removable Back and Arm, 2/cs
12025KD-1	Knocked Down, Tool-free Removable Back and Arm, 1/cs

SPECIFICATIONS

Description

Outside Legs	32.5" (W) x 24" (D)
Seat Dimensions	26.5" (W) x 16.5" (D)
Seat Height	18.5" - 23.5"
Weight	14 lbs.
Carton Shipping Weight (-1 -2)	14 lbs. 27 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

12011SS	1pr/bx
12011SL	1pr/bx

PARTS

Description

12025KDB	Back Assembly
12025KDBC	Back Clip
12025KDA	Arm
12025KDLEGS	Set of 4
12021C	Clip, 1/ea
12011SS	Sm Suction Cups Pr., 1pr/bx
12011SL	Lg Suction Cups Pr., 1pr/bx
12025KD9LB	Back Assembly
12025KD9LS	Seat Pad
12025KDC	Back Clip, 1/ea