



ITEM NUMBER: 081705375

ITEM NAME: Homecraft Bathtub Transfer Bench

MANUFACTURER PART NUMBER: GT-5006

OVERVIEW:

The Homecraft Transfer Bath Bench helps with safe transfers into the bathtub from a wheelchair. The bench offers a sturdy seat for the elderly or those with disabilities. This makes it easier to wash up without the risk of slips and falls.

The bench is made of a 1" sturdy aluminum frame that is corrosion resistant. The back can be reversed without tools for use on either side of the bath bench. Extra-large suction cups provide added safety during use. The bath bench has a weight capacity of 500 lbs.

Enjoy an easier shower with this bench. Prevent falls on slippery surfaces by sitting and sliding instead of stepping over the tub's side. It is also helpful for elderly parents or grandparents who lack the mobility to get into the bath or the strength to stand throughout the process.

SPECIFICATIONS:

Dimensions	32"W x 23"D x 18"-22"H
Brand	Homecraft
Weight	13.23 lbs
Material	Rust resistant aluminum frame and hard plastic Seat
Color	White
Latex Free	Yes
Base	32"W x 23"D
Capacity	500 lbs.
Floor to Seat	18"-22"H
Seat	26 1/2" W - 16" D
Seat to arm	9 1/2" H