



ITEM NUMBER: 7100182

ITEM NAME: Days Standard Aluminum Crutches - Youth

MANUFACTURER PART NUMBER: GT-1012YJMC2300HD

OVERVIEW:

Lightweight crutches feature push-button adjustment to provide a custom, comfortable fit for every user. Designed for home or clinical use, these crutches are ideal for keeping weight off an injured or post-operative foot, ankle, knee or leg. The crutches are sold in pairs and are latex free.

SPECIFICATIONS:

Material	Aluminum Crutches
Weight	3 lbs.
Weight Capacity	300 lbs.
Height	37"-46"
Warranty	Standard 1 year warranty on manufacturing defects
Cleaning	Can be cleaned with non-corrosive detergents or soap and water