

Size 1
Additional Information Leg Support: 10" (260 mm)
Back Support: 14"H (350 mm)
Head Support 6" (150 mm)
Overall Length 39" (975 mm)
Seat Depth 9" (225 mm)

Size 2
Additional Information Leg Support: 10" (260 mm)
Back Support: 16"H (400 mm)
Head Support 6" (150 mm)
Overall Length 42" (1050 mm)
Seat Depth 11" (275 mm)

Size 3
Additional Information Leg Support: 13" (340 mm)
Back Support: 18"H (460 mm)
Head Support 8.5" (215 mm)
Overall Length 52.5" (1312.5 mm)
Seat Depth 13" (335 mm)

Size 4
Additional Information Leg Support: 13" (340 mm)
Back Support: 22"H (550 mm)
Head Support 8.5" (215 mm)
Overall Length 59.5" (1487.5 mm)
Seat Depth 16" (400 mm)