

BV
MEDICAL

20-786-014-L
&
20-776-014-L

BV
MEDICAL

Clock Aneroid Sphygmomanometer Wall Mount, Adjustable, Black Nylon Cuff, Adult, Latex Free

Large, easy-to-read, 6" faceplate with white luminescent numerals
Adult size cuff, fits arm circumference range 23-40 cm (10 3/4" to 15 3/4") • Square plastic casing encloses a no stop pin, 300 mmHg gauge • Calibration tool allows for on-site calibration • LIFETIME calibration warranty +/- 3mmHg • 8' Latex Free coiled tubing (black) • Deluxe air release valve • Easy viewing with swivel feature

20-776-014-L



20-786-014-L

**Clock Aneroid Sphygmomanometer,
Rolling Stand Mounted
Adjustable Black Nylon Cuff
Adult
Latex Free
Adjustable floor stand (42.5" - 55")
Easy assembly Stand
Can be purchased separately**

BV Medical
28W206 Commercial Ave., Unit B
Lake Barrington, IL 60010
(888) 822-8293 | info@bvmedical.com
www.bvmedical.com