

# T3xm Treadmill



- Class IIa Med CE certified for therapeutic, rehabilitation and active-aging settings
- Intuitive LED console displays only the feedback essential to therapists and clinicians
- WiFi connectivity accommodates optional Matrix Asset Management system and Workout Tracking Network
- Integrated IR sensor stops belt motion when triggered by an object or person
- Low 0.16 km/h / 0.1 mph start speed in forward and reverse lets users with limited abilities take part in functional movement exercises
- Extended handrails with neutral-position, angled grips ensure secure entry and exit
- Extra-wide deck offers an added sense of security and freedom
- Maximum reverse speed of 6.4 kph / 4 mph
- Ultimate Deck System for superior performance and durability
- Matrix 4.2 HP AC Dynamic Response Drive System fine-tunes response based on footfall patterns for a smooth, consistent, natural workout

## CONSOLE

Display Type	Dot-matrix LED with profile display
Display Feedback	Time of Day, Time Elapsed, Time Remaining, Distance (Kilometers or Miles), Calories, Calories per Hour, Speed, Incline, Pace, Average Pace, Heart Rate, METs, Watts, Workout Profile, Forward/Reverse
Secondary Data Display	No
User-defined Multi-language Display	Accessible only through manager menu: English, German, French, Italian, Spanish, Dutch, Portuguese, Japanese, Swedish, Finnish, Turkish, Polish, Russian, Danish
Workouts	Manual, Rolling Hills, Fat Burn, 5k, 10k, Target HR, Gerkin Protocol, Army PFT, Navy PRT, Marine PFT, Air Force PRT, Physical Efficiency Battery (PEB), WFI Protocol
CSAFE Ready-FitLinxx™ Certified	Yes
IPTV Compatible	No
Pro:Idiom Compatible	No
FitTouch™ Technology	No
One-button Start	Yes
On-the-fly Program Change	Yes
Integrated Vista Clear™ Television Technology	No
FITCONNEXION™ Ready	No
WiFi Enabled	Optional
Bluetooth Enabled	No
iPod® / iPhone® Compatible	Charging only
Multimedia Playback	No
Personal Fan	No
Personal Trainer Portal Compatible	No
USB Port	Yes
Virtual Active™ Compatible	No

Asset Management Compatible	Yes
Workout Tracking Network Compatible	No
Web Connectivity	No
Facility Communication Portal/Calendar	No
RFID Compatible	No
Pause Function	Yes
Reading Rack	No

#### DRIVE SYSTEM

Drive System	Matrix 4.2 hp AC Dynamic Response Drive System with .1mA current leakage
--------------	--

#### FEATURES

Deck Type	Ultimate hard-wax reversible 1" deck
Belt Type	Habasit, 2-ply commercial grade
Running Area	152 x 51 cm / 60" x 20"
Step-on Height	19 cm / 7.5"
Cushion System	Ultimate Deck™ cushioning system
Incline Range	0-15% (590-kg / 1,300-lb. thrust elevation motor)
Speed Range	0.16 - 19.3 km/h / 0.1 - 12 mph forward / 0.16 - 6.4 km/h 0.1 - 4 mph reverse
Contact & Telemetric HR	Yes
Crossbar Controls	No
Tread Sense	Yes. Engages above 0.8 kph / 0.5 mph

#### NOTE

Speed Rating	24.1 km/h / 15 mph rating good for users up to 159 kg / 350 lbs.
--------------	--

#### TECH SPECS

Power Requirements	20-amp dedicated circuit required, non-looped grounded
Electrical Receptacle & Plug	NEMA 5-20R 120v
Assembled Dimensions	215 x 85 x 142 cm / 84.6" x 33.5" x 56"
Max User Weight	182 kg / 400 lbs.
Assembled Weight	172 kg / 378 lbs.
Shipping Weight	186 kg / 409 lbs.