

T655MD TREADMILL

The T655MD Series treadmill is a new addition to our medical treadmill line. The T655MD offers a greater incline range and long medical handrails that can rise, lower and angle inward. Additionally, a detachable step has been added, lowering the step-up height. The T655MD comes standard with a wired controller, allowing the physician or rehabilitation technician to easily control both speed and incline of the unit.

KEY FEATURES

- 5 HP AC-Servo motor
- ECO-Glide™ automatic lubrication system
- Detachable step
- -3% to 22% incline
- Reverse: 0.1-3 mph / Forward: 0.1-12 mph



TECHNICAL DETAILS

Unit Weight	464.2 lbs / 211 kg (with step)
Dimensions (LxWxH)	104 x 40 x 59.2 in / 264 x 101.6 x 150.4 cm (with step)
Running Area (LxW)	62 x 22 in / 158 x 56 cm
Speed Range	-3 to 12 mph / -4.8 to 19.3 kph
Incline Range	-3 to 22%
Drive Motor	5.0 HP
Power Requirements	15A / 120V / 60Hz 10A / 220V / 60Hz
Plug Requirements	NEMA 5-15P
Max User Weight	500 lbs / 227 kg
Step Up Height	With step: 6 in / 15.25 cm Without step: 10 in / 25.2 cm
Micro Amp Leakage Components	≤ 300µA
Readouts	Calories, Speed, Time, Distance, Cal/Hr, METs, Pace, Incline, Heart Rate, Weight Loss Zone, Cardio Zone, ActivZone
Workout Programs	Manual, Random, Interval, Glute, Hill, Track, Fat Burn, HRC Program, Bruce, Modified Bruce, Naughton

WARRANTY

Type	Full Commercial - no usage limit
Frame	Lifetime
Parts	5 Years
Wear Items	5 Years
Labor	3 Years

TECHNICAL DETAILS

Features	CSAFE port ECO-Glide auto lubrication Emergency stop Dual headphone jacks Detachable step Wireless heart rate receiver Adjustable height and angle handrails -3 to 22% incline Standard wired controller remote Integrated 3-speed fan Meets UL requirements for current leakage in medical settings
Optional Features	EDAN ECG Integration

CONSOLE



Tri-Color LED

SPECIAL FEATURES



Detachable step

SPECIAL FEATURES



Adjustable height and angle handrails