


KABOOTI®




New
Large Size
Available

MAKES
THE ORDINARY
DONUT
OBSOLETE


3-IN-1 COMPLETE COMFORT SEATING SOLUTION



COCCYX
CUSHION



SEATING
WEDGE



DONUT
RING

COCCYX CUTOUT



COMFORTS &
REDUCES
PRESSURE ON
YOUR TAILBONE

- ✓ Generous Coccyx Cutout
- ✓ Contoured for Comfort
- ✓ Tapered Front Edge
- ✓ Reduces Lower Back Pain
- ✓ Provides Stabilization
- ✓ Available in Blue & Gray



CONTOUR®

We Support You for Life

www.contourproducts.com

KABOOTI®



SCULPTED SURFACE
Comfortably fits your bottom's shape.

GENEREROUS COCCYX CUTOUT
Eliminates direct contact between you and the chair to relieve tailbone pressure during prolonged seating.

ERGONOMIC WEDGE DESIGN
Encourages proper posture.

MALE/FEMALE RECESS
Provides genital contact relief.


MACHINE-WASHABLE, STRETCH-KNIT COVER with carry handle
Allows for discreet use anywhere.

LARGE CENTER OPENING
Unique center opening removes pressure and connects to the coccyx cutout for cooling airflow.

INNOVATIVE STABILITY LOBES
Evenly distributes weight to provide greater stability and prevents tipping.

TAPERED FRONT EDGE
Reduces pressure points under thighs and allows legs to fit underneath desks, tables, and vehicle steering wheels.

The **KABOOTI®** is uniquely engineered to be superior to conventional donut cushions.



Pressure mapping* shows a comfortable seating with minimal pressure on the Kabooti, but the donut cushion shows extremely high pressure concentration (red and orange areas) on buttocks area. This can result in uncomfortable seating, reduced skin and surface blood circulation and potential damage to skin.

*Pressure maps performed with 175 lb. male, using XSensor pressure mapping equipment & software.

Product	Part #	UPC #	Case Dimensions	Case Pack
Kabooti - Blue Cover	30-751RB	737709003611	17"x 3" x 13.5"	4
Kabooti - Gray Cover	30-751RGY	737709003628	17"x 3" x 13.5"	4
Kabooti - Large - Blue Cover	30-740RB	737709004472	21"x 3" x 13.5"	4
Kabooti - Large - Gray Cover	30-740GRY	737709004489	21"x 3" x 13.5"	4

Email: sales@contourhealth.com

Contour Products, Inc.

Charlotte, NC

www.contourproducts.com | 1-855-896-8679

Invented and designed in the United States. Made in China.

Contour and Kabooti are registered Trademark used by Contour Products, Inc., under license from Banyan Licensing, L.L.C. US Patents: D690.537, D681.995. Other US and Foreign patents pending.

K_006_0718

© 2018. All rights reserved.