



A healthier,  
safer and  
more  
dignified  
walking  
experience.

Upright Posture Walker designed to be safer & more comfortable than traditional walkers and rollators.



Comfortable, integrated seat.

UPWalker is proudly produced by LifeWalker, a leading innovator in upright mobility.

**LifeWalker**  
Mobility Products

The revolutionary UPWalker® restores activity and independence for the elderly and mobility challenged.

- Enables users to stand upright and look ahead with better support, greater dignity and confidence to walk farther and remain more active.
- Preliminary clinical studies revealed improved walking posture, safer, more stable feeling, and reduced pain in legs, back, hands and wrists.
- Helpful in treatment and mobility for patients with Orthopedic, Cardiovascular, Neurological, Pulmonary and other health disorders.
- Stylish, easy-folding design, lightweight for transporting and storage, simple armrest adjustment for the perfect height in seconds.
- Cushioned sit-to-stand handles ease transition to standing position.
- Standard accessories include drink holder, seat backrest, personal bag.
- Sturdy frame with user weight capacity up to 350 lbs (Large model).

\* Limited Lifetime Warranty covers frame against defects in materials and workmanship.



**UPWalker®**

9530 Chesapeake Dr. Ste 501, San Diego, CA 92123 • (877) 488-0822 • LifeWalkerMobility.com  
UPWalker (patent pending) & LifeWalker are trademarks of ProtoStar, Inc., dba LifeWalker Mobility Products.



# UPWalker® features provide benefits not available in traditional walkers.



### Adjustable-Height Armrests

Padded armrests enable users of varying height to stand tall and look ahead with confidence and dignity.



### Lockable Brakes

Ergonomic handbrakes are easy to use for better control, stopping and parking.



### Comfortable Seat

Fabric seat with backrest support slides front to back for sitting or walking.



### Sit-To-Stand Handles

Padded handles make transition between sitting and standing easier and safer.



### Multi-Terrain Wheels

Large 8-inch rubber wheels help to move smoothly over bumps and uneven surfaces.



### Sturdy, Foldable Frame

Stylish, lightweight, folds easily for transport and storage yet strong enough to support 300 lbs.

## Helpful UPWalker Accessories:

- Backrest Support**  
 Detachable backrest provides comfortable support while sitting.
- Beverage Holder**  
 Attaches to handgrip arm for easy access to your water bottle.
- Personal Item Bag**  
 Polyester bag with fold-over top to conveniently stow personal items.
- Luxury Item Bag**  
 Elegant design, sturdy material, side pockets and easy-access zippered top.
- Shopping Bag**  
 Perfect for groceries, features handles to use as a carry bag.
- Smartphone Holder**  
 Adjustable to securely hold phones of various sizes for handsfree operation.
- Cane/Umbrella Holder**  
 Keep a cane or umbrella within easy reach with this two-piece holder.
- Flashlight/Taillight**  
 Serves as a headlight, flashing safety light and rear-facing taillight.

## 3 UPWalker Size Options for Best Fit

The UPWalker is available in 3 sizes to accommodate user height from 4'4" to 6'4" and weight up to 350 lbs.

### UPWalker Product Specifications:

FEATURES	SMALL (H200-S)	STANDARD (H200)	LARGE (H200-L)
<b>Product</b>			
User Height Range	4'4"-5'5" (132-165 cm)	4'7"-5'10" (140-178 cm)	5'0"-6'4" (152-193 cm)
User Weight Capacity	300 lb (136 kg)	300 lb (136 kg)	350 lb (159 kg)
Length	30.5" (77 cm)	33.5" (85 cm)	36" (91.4 cm)
Width	21.5" (54 cm)	23.5" (60 cm)	25.5" (65 cm)
Armrest Height Range	31.5"-40.5" (80-103 cm)	33"-45" (84-114 cm)	36"-49.5" (91-126 cm)
Seat Width	16" (40.5 cm)	18" (43 cm)	20" (51 cm)
Seat Height From Ground	20" (51 cm)	20" (51 cm)	21" (53.5 cm)
Sit-to-Stand Handle Height	28" (71 cm)	30" (76 cm)	32" (81 cm)
Unit Weight	21 lb (9.5 kg)	23 lb (10.5 kg)	25.5 lb (11.5 kg)
Wheel Size	8"		
<b>Folded Dimensions</b>			
Length	30.5" (77 cm)	33.5" (85 cm)	36" (91.5 cm)
Height	35" (89 cm)	37.5" (95 cm)	40" (101.5 cm)
Width	10.5" (27 cm)		
<b>Warranty</b>	Lifetime Warranty on frame; 6-months on non-durable components & accessories		

# UPWalker®

**SPECIAL OFFER!** Get 3 FREE Accessories and FREE Shipping!

Call (877) 488-0822 or visit LifeWalkerMobility.com to order an UPWalker, backed by a 30-Day Satisfaction Guarantee!