

Features and Benefits of the PATHWAY® CTS Pelvic Muscle Rehabilitation System

Features:	Benefits:
1. Documentation	1. Chart note and narrative reports, review and compare sessions, transfer to a specific location, PDF and Microsoft Word editable, compare sessions, and progress tracking
2. EMG, biofeedback and stimulation on one platform	2. Provide two channels of EMG and stimulation, monitor muscle contractions, and perform stimulation synchronously
3. Pre-set and modify protocols	3. Modify preset protocols enables you to store your protocol list with such things as: Set goals, change display type, single or split-screen, select colors, and designate markers
4. Premier privacy and compliance	4. Password protection, multi-level access, unique patient identifiers, and audit management included
5. Unlimited technical support with expanded hours	5. Phone and remote access, problem solvers, rapid response time
6. Accessories included	6. Starter kits include multiple accessories to begin use immediately
7. The PATHWAY® CTS-1500B includes laptop computer and printer	7. Customized for clinical use including pre-installed software for clinical use
8. File transfer	8. Automatic and manual database backups, database encryption, and file paths for document transfer
9. Treatment	9. Up to six treatments per patient session
10. Time tested products	10. Guaranteed quality, Made in the USA