

PATHWAY[®] CTS Pelvic Muscle Rehabilitation System Description

PATHWAY[®] CTS 1500 Pelvic Muscle Rehabilitation System, (Part # CTS-1500B)

Includes one Notebook Computer, One Color Printer, One USB Serial Interface Cable, One Pathway[®] CTS1500 Pelvic Muscle Rehabilitation System USB Flash Drive Software Program, One Pathway[®] CTS 1500 Pelvic Muscle Rehabilitation System Operator's Guide One Pathway[®] CTS 1500 Pelvic Muscle Rehabilitation System Module providing one Channel of Vaginal and Rectal Intracavity EMG or External Pelvic Muscle EMG, One Channel of Accessory Muscle EMG, One Channel of Vaginal and Rectal Electrical Stimulation

Starter Accessory Package

- (4) – Disposable Lead Wire Electrodes Sample Packets, Square Electrode, 24" Lead Wire, .060 Female Pin (7400)
- (4) - Pathway[®] Electrodes Sample Packets, (Rectangular, 3 Snaps in a Row Style), Use with Pathway[®] EMG Systems (6750)
- (1) – Ch A 10' Interface Cable for 7400 Disposable Lead Wire Electrodes (7300)
- (1) – Ch B Extended Pathway[®] Preamplifier, White, 10', Use with Pathway CTS 2000 or Pathway CTS 1500 (2583E)
- (4) - Pathway[®] Vaginal EMG/Stimulation Sensor (6330)
- (2) - Pathway[®] Rectal EMG/Stimulation Sensor (6340)

Support Services:

- (1) Unlimited hardware and software telephone Technical Support.
 - (2) Team Viewer
-

PATHWAY[®] CTS 1500 Pelvic Muscle Rehabilitation System, (Part # CTS-1500)

The PATHWAY[®] CTS 1500 **does not** include a **notebook computer** or **color printer**.

Minimum Computer Requirements

- Windows 10 Professional
 - 64-Bit Operating System
 - Intel i5 Processor or Higher CPU
 - RAM: 8 GB RAM
 - 500 GB Minimum Available Hard Drive Space (To Allow for Database Expansion)
-

The Prometheus Group[®]

One Washington Street, Suite 3171 Dover, NH 03820
Phone: 1 (800) 442-2325 | Fax: 1 (603) 749-0511 | theprogrp.com
Questions? Email us at performance@theprogrp.com

Share & Follow Us Online!



MKT-029 Rev. A.