

EMYO® 200B (Part #EM200B-N) Pelvic Floor & Orthopedic Rehabilitation System for Adult and Pediatric Physical Therapy

The EMYO® 200B system (Part #EM200B-N) contains the following items:

- (1) Notebook Computer
- (1) Wireless Mouse
- (1) Color Printer
- (1) 24" HDMI Monitor with Stand
- (1) USB Serial Interface Cable
- (1) EMYO® 200B USB Flash Drive Software Program includes Two Channel EMG Biofeedback with Games
- (1) EMYO® 200B Operator's Guide
- (1) EMYO® 200B Module providing two Channels of EMG
- (1) Carry Pouch
- (1) Heavy-Duty 9V Alkaline Battery

Starter Accessory Package

- (2) – Pathway® Preamplifier, Available in White or Gray, 5' (Part #2583)
- (1) – Pathway® Electrodes Sample Packet, (Rectangular, 3 Snaps in a Row Style) (Part #6750)
- (1) – Pathway® Adapter (Part #3660)
- (1) – Adapter for 7400 Disposable Lead Wire Electrodes (Part #7100)
- (4) – Disposable Lead Wire Electrodes Sample Packets, Square Electrode, 24" Lead Wire, .060 Female Pin (Part #7400)

Support Services:

- (1) Unlimited hardware and software telephone Technical Support.
 - (2) Team Viewer
-

The EMYO® 200B system (Part #EM200B) contains the following items:

The EMYO® 200B system (Part #EM200B) **does not** include a **notebook computer, color printer, or 24" extended monitor with stand**.

Minimum Computer Requirements

- Windows 10 Professional
 - 64-Bit Operating System
 - Intel i5 Processor or Higher CPU
 - RAM: 16 GB RAM
 - 500 GB Minimum Available Hard Drive Space (To Allow for Database Expansion)
-

The Prometheus Group®

One Washington Street, Suite 3171 Dover, NH 03820
Phone: 1 (800) 442-2325 | Fax: 1 (603) 749-0511 | theprogrp.com
Questions? Email us at performance@theprogrp.com

Share & Follow Us Online!



MKT-031 Rev. A.