



People.Health.Care.



Transelast® Classic

Soft, close-weave,
conforming gauze bandage.

Field of Application

- Compression bandage for lymphedema and edema
- Intended for bandaging of fingers, toes, and genitals
- Dressing retention

Features/Benefits

- Soft, comfortable, close-weave gauze bandage
- Ideal narrow width for bandaging fingers and toes
- Conforms well to body contours
- Air-permeable
- Extensibility is approximately 130%

Product composition

59% polyamide, 41% viscose

Not made with natural rubber latex



Transelast® Classic conforming bandage
rolled, bulk packed in a cardboard box

Size	Length		PH Ordering Info	
	unstretched	stretched	Item No.	Qty per box
4 cm / 1.6 in	1.7 m / 1.9 yd	4 m / 4.4 yd	7103175	20
6 cm / 2.4 in			7103176	20

