



118 Log Canoe Circle, Stevensville, MD 21666

Phone: 410-604-6393

Fax: 410-604-6398

www.3pointproducts.com

Product Improvement Notice

3pp® Carpal Lift™

(formerly 3pp Carpal Lift NP)

The 3-Point Products' 3pp Carpal Lift brace has been improved with the addition of stretchable binding around the body of the brace.

Binding offers a better fit, easier application and improved fabric stability for long lasting wear. A simple upgrade but one we are certain improves this functional support for those experiencing wrist pain and dysfunction due to ulnar or mid carpal instability and/or TFCC injury.

The 3pp Carpal Lift features:

- A unique strap and pad system applies counterforce support to relieve wrist pain
- Allows full finger and wrist motion for easy-to- wear for all day function
- Washable, breathable and latex-free foam-lined material
- Two sizes fit a wide range of hand sizes and is easily customized. Available in right and left hand models.
- Measure around the MP joints to determine size

Order Information

Size	Right	Left	Fits
Sm/Med	P2012-R23	P2012-L23	6" (15.2 cm) to 7.25" (18.4 cm)
Med/Lrg	P2012-R34	P2012-L34	7.25" (18.4 cm) to 9" (22.8 cm)

Improved product began shipping May 1, 2021

For access to updated images and product copy and to place orders, please contact

3-Point Products, Inc.

Service@3pointproducts.com

888-378-7763 M-F 8:30-4:30 EST