

# CHUDO



A Kismet Brand

[chudous.com](http://chudous.com)

## HELP YOUR PATIENTS RECOVER FASTER WITH HIVAMAT<sup>®</sup> DEEP OSCILLATION<sup>®</sup>

### Chudo Hivamat<sup>®</sup> Deep Oscillation<sup>®</sup>

Sports teams and therapists all over the world use Deep Oscillation<sup>®</sup>

Relieves pain, relaxes muscles, breaks up solids in the interstitium, and works on the lymphatic system to facilitate healing from all types of muscular micro trauma injuries. Learn at [www.chudous.com](http://www.chudous.com).

### The Authentic Next Generation Hivamat<sup>®</sup>



**Patented Technology**

US Patent# 7,343,203




**FDA Approved**

Registration Number: 3002892330



**Supported by decades of research and clinical application**



### Non-Invasive & Deep-Acting Clinically Proven Effects:

- Edema Reducing
- Reduce Swelling
- Pain Reducing
- Muscle Relaxing
- Anti-Inflammatory
- Decreases Inflammation
- Promotes Healing
- Improves Flexibility

**Electrostatic Oscillations Permeate Up To 8 cm Deep**



DOP = Deep Oscillation Personal

DOE = Deep Oscillation Evident

! = New and Improved


	<b>Richmar® Hivamat® (discontinued)</b>	<b>Chudo Hivamat® (next generation)</b>
MODEL	Hivamat® 200 Portable Hivamat® 200 Evident	Chudo Personal (DOP) #02040P-USA-CHUDO Chudo Evident (DOE) #01130-USA-CHUDO
FREQUENCY RANGE	Portable: 5 - 250 Hz Evident: 5 - 250 Hz	Personal: 5 - 250 Hz Evident: 5 - 250 Hz
CHANNELS	Portable: One channel Evident: Two channels	Personal: One Channel Evident: two Channel
DISPLAY	Portable: Monochromatic Evident: Multicolor	Personal: Multicolor ! Evident: Multicolor
POWER	Portable: 1.2V NiMH RCHRGBL Battery (4X) Evident: 100 - 240 VAC	Personal: 7.2V Li-Ion 24 Wh RCHRGBL Battery (1x) ! Evident: 100 - 240 VAC
THERAPY CARD	Three: separate practice-specific treatment protocol cards	One: one card includes all treatment protocols !
FREE PARAMETER SETTINGS	Portable: No Evident: Yes	Personal: Yes ! Evident: Yes
ANATOMICAL ANIMATIONS	Portable: No Evident: Yes	Personal: Yes ! Evident: Yes
DESCRIPTIVE TREATMENT NAMES	Portable: No Evident: No	Personal: Yes ! Evident: Yes !
SERVICE QR CODE	Portable: No Evident: No	Personal: Yes ! Evident: Yes !
FDA-APPROVED	Yes	Yes



**In-Clinic Professional Use  
Chudo Evident**



**In-Clinic & At-Home Use  
For Professionals and Patients  
Chudo Personal**



**Chudo Hivamat Deep Oscillation by Kismet**