

REBOUND® POST-OP ELBOW

The Rebound Post-Op Elbow brace by Össur combines intuitive features for fitters and users with a design that aims to improve the post-op experience. Combining a drop-lock hinge and four telescoping paddles, individual post-op protocols are easily catered to. The Rebound Post-Op Elbow, lightest of its kind, promotes user comfort with straps that stay in place, anatomical left/right fit with moldable stability cuffs, and blue patient touch-points for easy education.



INDICATIONS FOR USE

Indications requiring elbow immobilization and/or range-of-motion restriction. This may include:

- Stable fractures of the elbow, distal humerus, proximal radius or ulna
- Tendon and ligament injuries or repairs
- Collateral ligament reconstructions
- Elbow hyperextension
- Triceps tendon repair
- Chronic elbow injuries
- Tennis elbow release

FEATURES AND BENEFITS

- Designed to be the most lightweight brace on the market.
- Telescoping paddles relocate to individualize length and keep incision sites clear.
- Easy education for users with intuitive patient touch-points, including Blue is You™ buckles and numbered strapping system.
- Drop-Lock hinge enables secure range-of-motion with bump-resistant pull adjustment stops and optional ROM lock-outs.
- Moldable stability cuffs designed to be individually customized for improved fit and support.
- Strap guides maintain strap location for added security and comfort.
- Universal size, left and right version.

REQUIRED MEASUREMENTS

- Universal, one size

REBOUND® POST-OP ELBOW

ORDERING INFORMATION

Side	Part#
Left	B-8391680
Right	B-8391580

ACCESSORY KITS

Description	Part#
Elbow Wrist Control Kit	B-8390003
Elbow ROM Locks Kit	B-8390005

REPLACEMENT PART KITS

Description	Part#
Elbow Liner Kit	B-8390002
Elbow Sling Kit	B-8390004
Strap Tab Kit	B-8390006
Adhesive Hook Sticker Kit	B-8390007
Anti-Migration Pad Kit	B-8390008
Elbow Strap Kit	B-8390009