

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



Designed to support the busy practitioner and multiple patients, solid hardwood construction means our Classic Wood Mat Table is built for strength and stability during practitioner-assisted exercise and other treatments. Our 4 Leg Mat table comes standard with solid hardwood on the base frame, four legs and supports with a thick plywood top construction. The mat table also features a layer of comfortable 1.5" high density foam with softly rounded corners and triple-bolted leg assembly for extra durability.

STANDARD FEATURES:

- Load Capacity: 450 lbs.
- Solid hardwood frame and 4 legs
- 3/4" Plywood top
- 1.5" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- ADA Compliant

Classic Wood
MAT TABLES



MAT-4872-4L

STANDARD TABLES:

- | | | |
|--|--|--|
| <input type="checkbox"/> MAT-4872-4L
48"W x 72"L x 19"H
Qty _____ | <input type="checkbox"/> MAT-4884-4L
48"W x 84"L x 19"H
Qty _____ | <input type="checkbox"/> MAT-6084-4L
60"W x 84"L x 19"H
Qty _____ |
|--|--|--|

STANDARD WOOD CHOICES:

- | | |
|--|--|
| <input type="checkbox"/>  Oak | <input type="checkbox"/>  Maple |
|--|--|

WOOD STAIN OPTION:

- MAT-Stain

STAIN CHOICES:


- | | |
|--|---|
| <input type="checkbox"/>  Warm Copper | <input type="checkbox"/>  Rustic Fall |
| <input type="checkbox"/>  Deep Garnet | <input type="checkbox"/>  Warm Caramel |
| <input type="checkbox"/>  Fossil Gray | <input type="checkbox"/>  Black Beauty |
| <input type="checkbox"/>  Tungsten | <input type="checkbox"/> Custom Wood Stain Color
_____ |

OPTION:

- LOGO-BB** - Table Sling



STANDARD UPHOLSTERY CHOICES:

- | | |
|--|--|
| <input type="checkbox"/>  Royal Blue | <input type="checkbox"/>  American Beauty |
| <input type="checkbox"/>  Imperial Blue | <input type="checkbox"/>  Lark |
| <input type="checkbox"/>  Blue Ridge | <input type="checkbox"/>  Burgundy |
| <input type="checkbox"/>  Deep Violet | <input type="checkbox"/>  Dove |
| <input type="checkbox"/>  Forest | <input type="checkbox"/>  Charcoal |
| <input type="checkbox"/>  Black | <input type="checkbox"/>  Mocha |
| <input type="checkbox"/>  Mandarin Orange | <input type="checkbox"/>  Putty |