

AA2116 Without Lid 150mm (6")
AA2116L With Lid 150mm (6")

AA2114Y Without Lid 100mm (4")
AA2114L With Lid 100mm (4")

AA2112 Without Lid 50mm (2")
AA2112L With Lid 50mm (2")



UK contact details

Customer Services
Tel: 08444 124 330
Fax: 08448 730 100
Email: retail.sales@pattersonmedical.com

International contact details

Customer Services
Tel: +44 (0) 1623 448 706
Fax: +44 (0) 1623 448 784
Email: international@pattersonmedical.com

Visit us on-line at:

www.pattersonmedical.co.uk

Your local Savannah stockist is:

Blank box for local stockist information.



Patterson Medical Ltd.
Nunn Brook Road
Huthwaite
Sutton-in-Ashfield
Nottinghamshire
NG17 2HU, UK

2014 Model



Savanah
a Patterson Brand

Raised Toilet Seat

Savannah Raised Toilet Seat



Gently contoured surface - improves comfort and seating position



Front and rear cut outs
- for convenient personal cleaning



Clip-in lid system - easy to remove and attach

Back to front downward slope - assistance for movement on and off the toilet



Smooth aperture - no dirt traps and easy to clean

Spares & Optional Extras

09 154 3305	Pair of spare brackets
09 154 3313	Pair of spare seat hinges
09 154 3321	Pair of spare rear caps
09 154 3339	Optional lid with hinges
AA2602	Savannah Splash Sentry

Specifications:

Available in three height sizes, without or with Lid.

Overall length & width:

Length 407mm (16")

Width 370mm (14½") excluding brackets

Aperture:

Length 266mm (10½")

Width 207mm (8")

Maximum user weight:

190kg / 30st



Larger hand wheel - quicker and easier to adjust



Adjustable seat brackets - secures the seat in position