

# TENS 3000™

## Analog Unit, Three Mode

The TENS 3000 uses a microprocessor which gives exact dial performance and are as accurate as a digital unit, yet easy to operate. Benefits include maximum full strength power of 80 mA. The safety amplitude cap protects the amplitude knobs from accidental bumps or knob movements.

### FEATURES:

- System Includes Device, Lead Wires, Four Self-Adhesive Reusable Electrodes, 9V Battery, Hard Plastic Carrying Case and Instruction Manual
- One Year Warranty

### TECHNICAL SPECIFICATIONS:

- **Power Source:** 9-volt battery
- **Modes:** Normal, Modulate and Burst
- **Channels:** Dual
- **Timer:** 15 and 30 minutes or continuous
- **Pulse Width:** 50-250 Microseconds, Adjustable
- **Pulse Rate:** 2-150 Hz
- **Pulse Amplitude:** 0-80 mA each channel
- **Waveform:** Asymmetrical biphasic square pulse

### ORDER DETAILS:

7025239

Richmar TENS 3000  
Analog Unit, Three Mode

