

EMS 7500

Digital EMS

The EMS 7500 Digital Device features a patient compliance meter, three modes, safety amplitude cover, sturdy metal belt clip, rubber side railings for maximum handling and a premium silver finish.

FEATURES:

- System Includes Device, Lead Wires, Four Self-Adhesive Reusable Electrodes, 9-Volt Battery, Hard Plastic Carrying Case and Instruction Manual
- One Year Warranty
- RX Only

TECHNICAL SPECIFICATIONS:

- **Patient compliance meter** can store up to 60 sets of operation records and 999 hours
- **Power Source:** 9-volt battery
- **Channel:** Dual, isolated between channels
- **Pulse Width:** 50-300 μ s, adjustable
- **Pulse Rate:** 2-120 Hz, adjustable
- **Pulse Ramp:** 1-8 seconds
- **Pulse Amplitude:** 100 mA peak into 500 phm load
- **Contraction Time:** 0-90 seconds
- **Relaxation Time:** 2-90 seconds
- **Waveform:** Asymmetrical Bi-Phase Square Pulse
- **Timer:** 1-60 minutes or continuous



ORDER DETAILS:

DE7502

EMS 7500 - Digital EMS