


CLEANPatch[®]

Instantly restores damaged
mattress covers
to an intact and hygienic state



Save Time



Cost Effective



Sustainable

Why use CleanPatch®?

CleanPatch® repairs **hospital mattresses** and **stretchers**. Incorporating CleanPatch® and surface repair into your equipment maintenance program will save you time and money while reducing waste, all without compromising patient safety.

Clinically Tested

- Used by thousands of healthcare facilities across the globe
- Early intervention reduces the risk of cross contamination

Cost Effective

- Save up to 95% on equipment replacement
- Patented peel and stick technology repairs damage in seconds
- Eliminate administrative costs and time associated with equipment replacement


Sustainable

- Reduce landfill waste by preventing unnecessary equipment replacement

Product Information

Small

5.1 cm/
2 inches




Damage size:

up to 1.3 cm/0.5 inches

There must be a >1.3 cm/0.5 inch margin around the damage site to ensure optimal adhesion.

Medium

8.9 × 8.9 cm/
3.5 × 3.5 inches




Damage size:

up to 6.4 cm/2.5 inches

Large

8.9 × 15.2 cm/
3.5 × 6 inches









Damage size:

up to 12.7 cm/5 inches

Available in various packaging configurations.
Please contact us or your local distributors for more information.

CleanPatch® is an FDA, Health Canada, UKCA, and CE Registered Class 1 medical device.

Instructions



Surface Inspection and Repair Certification Program

- A virtual training program designed to incorporate surface repair into your equipment maintenance program

Contact us to implement a surface repair program at your facility:

1-888-623-7085 • info@surfacemedical.ca • cleanpatch.ca/training

CLEANPatch[®]

Visit cleanpatch.ca

for a list of distributors and approved service providers,
and to learn more about other **CleanPatch[®]** products.

Awards and recognition



ERNEST C. MANNING
INNOVATION
AWARDS



START UP
CANADA.



ASTech
Foundation
Inspiring Innovation

