

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



Comprised of 4 or 6 sections, the Thera-P Traction Tables use electric power to help with manual adjustment, allowing for ease of use and control of table positioning. For the practitioner, the Thera-P Traction Table offers a 20-37" full height adjustment and section mobility so the patient can be placed in the proper position for maximum results with minimum effort. Each section is adjustable, providing numerous achievable patient positioning options.

STANDARD FEATURES:

- Load Capacity: 650 lbs.
- 27"W x 78"L x 20-37"H
- Durable and stylish steel base
- White powder coat base
- Whisper quiet lift system
- Height adjustment: 20-37"
- Dual foot bar control
- 27" Wide cushions in 2 different top configuration options
- 2" High density foam cushions
- Rounded corners
- Friction-free gas shock cushion separation with lock and unlock lever
- (2) Low profile locking casters
- Head section angle range: +30° to -90°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- Lift back: 0 to 70°
- Table comes standard with traction head arm and adaptable mounting plate

Thera-P™
TRACTION TABLE



TXT-1004

TABLE CONFIGURATIONS:

TXT-1004


4 Section
Table Dimension: 27"W x 78"L x 20-37"H
Overall Dimension: 110"W x 64"L x 20-37"H

Qty _____

TXT-1006

6 Section
27" Width
Table Dimension: 27"W x 78"L x 20-37"H
Overall Dimension: 110"W x 64"L x 20-37"H

Qty _____



OPTIONS:

- TXT-ACC Accessories Kit**
Traction Head
Cervical Attachment
Traction Belt
Knee and Cervical Bolster Set
(1) Cervical Head Pillow
- S1354 - Strapping System**
Traction Belt
Knee and Cervical Bolster Set
(1) Cervical Head Pillow
- S1353 - Traction Belt**

TRACTION HEAD FEATURES:

- Crystal clear displays
- Static, Intermittent, and Cyclic Traction
- Variable speed control
- Minimum and maximum force settings 2.5kgs/5.5lbs - 91kgs/200lbs
- Patient safety switch with audible alarm
- Traction force compensation
- Audible alarm protection monitoring of traction forces
- Continuous electronic monitoring of traction forces



*Knee & Cervical
Bolster Set and
Traction Belt*



Traction Head



*Cervical
Attachment &
Head Pillow*

STANDARD BASE COLOR:

- White

CUSTOM BASE COLORS:

- | | |
|--|--|
| <input type="checkbox"/> <input type="checkbox"/> Black | <input type="checkbox"/> <input type="checkbox"/> Green |
| <input type="checkbox"/> <input type="checkbox"/> Light Gray | <input type="checkbox"/> <input type="checkbox"/> Forest Green |
| <input type="checkbox"/> <input type="checkbox"/> Charcoal | <input type="checkbox"/> <input type="checkbox"/> Purple |
| <input type="checkbox"/> <input type="checkbox"/> Royal Blue | <input type="checkbox"/> <input type="checkbox"/> Yellow |
| <input type="checkbox"/> <input type="checkbox"/> Navy Blue | <input type="checkbox"/> <input type="checkbox"/> Orange |
| <input type="checkbox"/> <input type="checkbox"/> Red | <input type="checkbox"/> <input type="checkbox"/> Brown |
| <input type="checkbox"/> <input type="checkbox"/> Burgundy | <input type="checkbox"/> <input type="checkbox"/> Bright Blue |

STANDARD UPHOLSTERY CHOICES:

- | | |
|---|---|
| <input type="checkbox"/> <input type="checkbox"/> Royal Blue | <input type="checkbox"/> <input type="checkbox"/> American Beauty |
| <input type="checkbox"/> <input type="checkbox"/> Imperial Blue | <input type="checkbox"/> <input type="checkbox"/> Lark |
| <input type="checkbox"/> <input type="checkbox"/> Blue Ridge | <input type="checkbox"/> <input type="checkbox"/> Burgundy |
| <input type="checkbox"/> <input type="checkbox"/> Deep Violet | <input type="checkbox"/> <input type="checkbox"/> Dove |
| <input type="checkbox"/> <input type="checkbox"/> Forest | <input type="checkbox"/> <input type="checkbox"/> Charcoal |
| <input type="checkbox"/> <input type="checkbox"/> Black | <input type="checkbox"/> <input type="checkbox"/> Mocha |
| <input type="checkbox"/> <input type="checkbox"/> Mandarin Orange | <input type="checkbox"/> <input type="checkbox"/> Putty |

CUSTOM BASE OPTIONS:

- S1178**
Table Base Graphics
- PT9020-CC**
Pick Custom Base Color
- S2000**
Custom Table Base Powder Coat Color and
Table Base Graphics

