

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



Comfortable, safe, and durable, the Elevating Therapy Tilt Table is ideal for nursing home and rehabilitation facility use. Smooth and secure elevation adjusts from 22.5-39.5" high, and tilts from a horizontal position to 90°. A heavy-duty powder coat frame supports a load capacity of 350 lbs. and includes an extra-large foldaway foot plate with non-skid surface and dual safety straps for added stability. The 2" high density foam cushions with two configurable options are designed for patient comfort. The two section cushion configuration includes a manual lift back that adjusts to 70°.

STANDARD FEATURES:

- Load Capacity: 350 lbs.
- 27"W x 82"L x 22.5-39.5"H
- White powder coat base
- (2) Safety straps
- Extra large foldaway foot plate with non-skid tape
- Hand control for tilt
- Single foot control for lift
- Tilt angle indicator: 0 to 90°
- 27"W Cushions with 2 different top configuration options
- 2" High density foam cushion
- Rounded corners
- Nose hole with plug
- (4) Individual locking casters
- Manual lift back: 0 to 70°
- Accommodates patient lift system

Elevating Therapy
TILT TABLE

HY1002



TABLE CONFIGURATIONS:

- | | |
|---|---|
| <input type="checkbox"/> HY1001
1 Section
27" Width
Qty _____ | <input type="checkbox"/> HY1002
2 Section
27" Width
Qty _____ |
|---|---|

STANDARD BASE COLORS:

- White

CUSTOM BASE COLORS:

- | | |
|-------------------------------------|---------------------------------------|
| <input type="checkbox"/> Black | <input type="checkbox"/> Green |
| <input type="checkbox"/> Light Gray | <input type="checkbox"/> Forest Green |
| <input type="checkbox"/> Charcoal | <input type="checkbox"/> Purple |
| <input type="checkbox"/> Royal Blue | <input type="checkbox"/> Yellow |
| <input type="checkbox"/> Navy Blue | <input type="checkbox"/> Orange |
| <input type="checkbox"/> Red | <input type="checkbox"/> Brown |
| <input type="checkbox"/> Burgundy | <input type="checkbox"/> Bright Blue |



HY1002

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|--|
| <input type="checkbox"/> Royal Blue | <input type="checkbox"/> American Beauty |
| <input type="checkbox"/> Imperial Blue | <input type="checkbox"/> Lark |
| <input type="checkbox"/> Blue Ridge | <input type="checkbox"/> Burgundy |
| <input type="checkbox"/> Deep Violet | <input type="checkbox"/> Dove |
| <input type="checkbox"/> Forest | <input type="checkbox"/> Charcoal |
| <input type="checkbox"/> Black | <input type="checkbox"/> Mocha |
| <input type="checkbox"/> Mandarin Orange | <input type="checkbox"/> Putty |