



# rolyan<sup>®</sup>

for Hand Therapy.

MADE IN THE USA



## Rolyan Splint Bath

Introducing a not-so-basic splint bath with innovative features and a modern look

- ✓ **Modern look**  
Made of a thermoset material known for its insulation, flammability resistance and durability. Features a stainless steel hinged bi-fold lid for easy access from three sides
- ✓ **Auto on/off**  
Automatically preheats the water to reach a set temperature on a set day and time. Ready when your clinic opens and energy efficient at night when not in use
- ✓ **Heats in 25 minutes**  
Internal heat exchanger is immersed directly in the water basin, ensuring fast heating and uniform heat distribution

### Key features

- Auto on/off
- Fast and uniform heating
- User sets time and temperature
- Audible/visual alerts
- Low water sensor
- Rapid water evacuation
- Bi-fold stainless steel lid
- Digital touch interface

Rolyan is a registered trademark of Performance Health and/or its subsidiaries and may be registered in the United States and other countries. © 2023 All rights reserved. Unauthorized use is strictly prohibited.

**7026425**  
**Rolyan Splint Bath**

Tub dimensions:  
12" X 19" X 6"



Handy storage



Bi-fold stainless steel lid



Rolyan, where ingenuity and artistry go hand in hand.