

360 Air Walker



Key Features:

- **Semi-rigid shell** – offers maximum support to the leg/ankle, ensuring maximum comfort
- **Adjustable toe guard** – provides additional protection for the foot
- **Shock absorber insole** – reduces the impact on the foot when walking
- **Inflatable air bladders** – increase comfort and reduce swelling
- **Rocker sole** – helps promote a natural gait

Size	Men's shoe size	Women's shoe size	Short item #	Tall item #
XS	–	3-5	7026976	7026971
S	5-6	5.5-7	7026977	7026972
M	6.5-10	7.5-11	7026978	7026973
L	10.5-12	11.5-13	7026979	7026974
XL	12.5-13	–	7026980	7026975

Learn More



Sammons Preston is a registered trademark of Performance Health and/or its subsidiaries and may be registered in the United States and other countries.
© 2024 All rights reserved. Unauthorized use is strictly prohibited. SPWB1123 (11/23)